





# Life Is On Schneider

## Schneider Electric ANALECT® RefinIR™

NextGen Hydrocarbon Analyzer for Lab Automation Process Development The all new ANALECT RefinIR is a fully integrated laboratory autosampler and FTIR instrument. It is designed to measure liquid hydrocarbons including gasoline, diesel and crude oils in a single autosampler. It features the ANALECT Diamond 20<sup>™</sup> Transept interferometer platform with the improved DCM 600 indicate acquisition electronics. This provides improved performance utilizing SpectraQ software with a customizable HMI for operating the system. Determination of physical properties and chemical composition can be achieved with the touch of a button. The RefinIR system automates calibration data collection in support of the Hydrocarbon SmartSystem® on-line process analyzer.

#### Spectra Q<sup>™</sup> Software Simplifies the Data Collection Process

With one touch of the button, SpectraQ software allows the user to automatically collect FTIR spectra of a full carousel of samples. In addition, the software provides the user the ability to perform a full ASTM D6122 validation to ensure spectral integrity.

- Full chemometric modeling capability including SpectraQuant is available.
   Optional third party chemometric software can be added including MATLAB, Eigenvector SoloPredictor and Unscrambler.
- Communications options Modbus, OPC and ethernet to LIMS.

# Spectra Q

### Specifications

- Optibus™optical transfer components integrate with Diamond 20™ Transept interferometer
- · Single-cell system for analyzing liquid hydrocarbons
- 0.25 and 0.5 mm pathlength liquid transmission cell
- Maximum vial-to-vial cross contamination : < 0.1%</li>
- · Thermo-electrically stabilized DTGS detector
- 115 VAC/60 Hz or 230 VAC 50/60 Hz
- Automated sampling system including carousel for twenty-four (24) 15ml vials
- Two (2) liter vessels for wash and validation solvents
- Dimensions: 74H x 60W x 45D cm ( 29H x 24W x 18D in.) Weight: 90 kg (200 lb)
- Oven controllers and intrinsically safe sensors to enable the sampling of flammable and non-flammable samples over the range of 40°C +/- 1°C to 125°C +/-1°C
- Analog input card for the host computer to read cell temperature

Schneider Electric Analyzer Systems (Formerly AIT) 2121 Aviation Drive Upland CA 91786 USA www.AITanalyzers.com





## RefinIR Features

- RefinIR delivers automated spectral data acquisition of multiple liquid re finery stream samples from crude oils to blended gasoline and diesel.
- Autosampler can be heated up to 125°C to allow sampling of heavy crudes and intermediates, such as FCC feed.
- Utilizes the ANALECT XCORR<sup>™</sup>
  model transference protocol allowing data to be transferred to the
  on-line ANALECT Hydrocarbon
  SmartSystem.<sup>®</sup>
- Holds 24 15 ml vials and automatically transfers the sample to a flowthrough transmission cell embedded in the autosampler.
- Includes compliance to ASTM
   D6122 via a wash and validation
   functionality that eliminates cross contamination of samples and pro vides a quick assurance of spectral
   quality.
- Reduces the laboratory benchtop foot print of the prior ANALECT RefinIR by up to 60%.
- Newly designed carousel for robust operation.